

EMCLAB Bottle Top Dispenser

Features:

- Quick, easy and precise volume adjustment
- High quality materials like PTFE, PFA, FEP, borosilicate glass, platinum-iridium with high chemical resistance again acids, solvents and bases. Please check before use!
- Positive displacement plunger – no tools necessary
- Easy to dismount for cleaning
- Completely autoclavable up to 121 °C
- Certified of conformity (DIN 12600)
- Margin of error to DIN EN ISO 8655-5: A < +/- 0.5%
- 3 Adapters for different GL-threads included

EMCLAB Bottle Top Dispenser

Article No.	Volume (ml)	Price €	Special Price €
EMC-20054010	0.5 - 2.5	160.00	120.00
EMC-20054020	1 – 5	160.00	120.00
EMC-20054030	2 – 10	160.00	120.00
EMC-20054040	5 – 25	178.00	132.00
EMC-20054050	10 – 50	186.00	138.00

Special offer

EMCLAB Bottle Top Dispenser

Article No.	Volume (ml)	Increment (ml)	Thread (GL)	Adapter (GL)
EMC-20054010	0.5 - 2.5	0.05	32	28, 40, 45
EMC-20054020	1 – 5	0.20	32	28, 40, 45
EMC-20054030	2 – 10	0.50	32	28, 40, 45
EMC-20054040	5 – 25	0.50	45	32, 38, S40
EMC-20054050	10 – 50	1.00	45	32, 38, S40

Prices: ex works Duisburg, Incoterms 2010

Special offer is valid as long as stocks last (without glass)

Minimum purchase order: 2 pieces

