

Introducing the Revolutionary Magnetic Quick Change Cell Holders


Upgrade your UV/VIS spectrophotometry workflow with our new Magnetic Quick Change Cell Holders. Using a unique magnetic system, this solution allows quick transitions between different cell sizes without tools or realignment.

Key benefits include:

- Fast & Flexible: Switch cell sizes effortlessly with no adjustments needed.
- Batch Convenience: Swap pre-prepared sets of 4 or 8 cuvettes at once.
- Perfect Alignment: Magnetic design ensures precise beam centering every time.
- Smart technology: Built in NFC chip for marking and identification.

Exclusively developed by us, this revolutionary holder is the first of its kind. Say goodbye to the hassle of traditional cell holders and step into a new era of efficiency and accuracy.

See the following technical drawing to enhance your workflow and increase efficiency today.



Cell Holder	Description	Art. No.
	Slider 4 - magnetic slider for easy installing Cuvette slider for 4-position Quick Change cell holder 10x10 mm, 10x50 mm, 10x100 mm, without Quick Change cell holder, incl. cell holder rod	EMC-00301
	Slider 8 - magnetic slider for easy installing Cuvette slider for 8-position Quick Change cell holder 10x10 mm, without Quick Change cell holder, incl. cell holder rod	EMC-00302
	4-Position Quick Change Cell Holder 10x10 For four cuvettes of the size 10x10 mm	EMC-00304
	8-Position Quick Change Cell Holder 10x10 For eight cuvettes of the size 10x10 mm	EMC-00305
	4-Position Quick Change Cell Holder 10x50 For four cuvettes of the size 10x50 mm	EMC-00306
	Cell Insert 20 Top cover for 4-position Quick Change cell holder 10x50 mm, for usage of 10x20 mm cuvettes	EMC-00303
	4-Position Quick Change Cell Holder 10x100 For four cuvettes of the size 10x100mm	EMC-00307

Including programmable NFC chip.
 Additional cell holders on request.